



THD SLIDE DOPPLER PROBE

800056	800063
NOT STERILE - REUSABLE	STERILE - DISPOSABLE
	
5 Shot (90 min/Shot)	One Shot – 90 min

Use

THD Slide Doppler Probes are intended to be used with THD Revolution device, for surgical treatment of haemorrhoidal disease, using the Transanal Haemorrhoidal Dearterialisation technique.

General description

THD Slide Doppler Probes shall be used with THD Slide kits and have been designed for maximum sensitivity to identify blood vessels.

Its small size allows the physician easy handling and when used with THD Slide anoscope/proctoscope, which has a port to lodge the Doppler, the operating field is completely unobstructed.

THD Slide Doppler Probes are available in both reusable (for limited treatments) and disposable versions.

THD Disposable Probe is individually packaged and STERILE.

THD Disposable Probe may be sterilised by Ethylene Oxide (EtO).



Dimensional Features	
THD Reusable Doppler – THD Disposable Doppler	
	
Body Dimension	14x70x7 mm
Cable Length	2000 mm

Information		
REF	TRADE NAME	Sterilisation - Resterilisation
800056	THD DOPPLER PROBE 5s/NS	ETO, Plasma Gas (Hydrogen-peroxide, using Sterrad system) or Peracetic acid using Steris system. * Do not sterilize in an autoclave or at temperatures exceeding 70°C.
800063	THD DOPPLER PROBE 1s/ST	Re-sterilization – NOT ADMITTED

* Sterilization method should be validated by final user.

THD Slide Doppler Probe			
REF	Qty	DESCRIPTION	EMDN
800056	1	THD DOPPLER PROBE 5s/NS – Not Sterile, 5 shots	Z12059005
800063	1	THD DOPPLER PROBE 1s/ST – Sterile, disposable, one shot	Z12059005