



# ROLLING CART FOR THD REVOLUTION



## Use

The Rolling Cart creates a mobile stand during procedure and storage of the THD Revolution (doppler dearterialiser unit).

## General Description

The Rolling Cart creates greater mobility, ease of orientation and negates the need to move other instruments in the room.

The cart is equipped with a handle for better control and ease of use.

Dimensional Features	
Dimension	Height: 1000 mm (39 in); Pivoting diameter 610 mm (24 In)
Weight	13 kg (28.7 lb)

THD ROLLING CART		
Qty	P/N	Description
1	880014	Rolling Cart for THD Revolution